

Model	Go Chair®
Part Number	GoChair
Weight Capacity	300 lbs.
Maximum Speed ¹	Up to 3.8 mph
Ground Clearance ³	1.75" at motor
Turning Radius ³	25.5"
Overall Length ³	28.5" (33" with footplate)
Overall Width ³	25.5"
Front Wheels	6" solid casters
Drive Wheels	9" solid
Rear Casters	2" solid anti-tips
Suspension	Limited
Range Per Charge ^{1,2}	Up to 11.2 miles
Heaviest Piece When Disassembled ³	36 lbs. (rear section)
Standard Seating	Type: Compact, foldable Material: Black vinyl Dimensions: Width: 18", Depth: 17", Weight: 30 lbs.; Width: 20", Depth: 18", Weight: 32 lbs.;
Dual Braking System	Regenerative and electromechanical
Standard Electronics	40 amp, PG GC3
Battery Requirements ⁴	20 Ah battery pack Weight: 31.4 lbs.
Battery Charger	Off-board, 2 amp
Warranty	Frame: 5-year limited Drivetrain: 13-month limited Electronics: 13-month limited Batteries: 13-month limited

HCPCS Code: K0899

**POWER MOBILITY DEVICE, NOT CODED BY DME PDAC
OR DOES NOT MEET CRITERIA**

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- Battery weight may vary based on battery manufacturer.

PRIDE® MOBILITY

800-800-8586 • pridemobility.com



1. FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

061826

GO CHAIR

LIGHTWEIGHT TRAVEL AND INDEPENDENCE ON THE GO



FDA Class II
Medical Device¹

800-800-8586 • pridemobility.com



GET OUT THERE

WITH THE #1 PORTABLE WHEEL CHAIR



The Go Chair® is #1 for portability! With convenient feather-touch disassembly, the Go Chair allows you to experience lightweight travel and independence on the go. Enjoy dual swing-away storage bins, an array of contemporary colors, and a compact turning radius for outstanding maneuverability.

FEATURES

- 18" x 17" or 20" x 18" standard 5-point swivel seat with latch release
- Compact turning radius
- Dual swing-away, under seat storage bins
- Height and width adjustable armrests
- One hand feather-touch disassembly
- 60" lap belt
- The heaviest piece is only 36 lbs.

OPTIONS

- Medium back memory seat options available in 18" x 18" or 20" x 18"
- Inline right and left swing-away joystick

AVAILABLE COLORS



OPTIONAL ACCESSORIES*

- Cane/crutch holder (3" single, 3" double)
- Cup holder
- Cell phone holder
- Lap belt (50", 70")
- Oxygen tank holder
- Rear basket
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger

MILITARY PATCHES



* Only one rear accessory can be used at a time.



GC3 controller with quick disconnect.

NEW rear basket for convenient storage

Optional swing-away joystick swings out at an angle for convenient access