

Jazzy® EVO 614HD¹	Group 2: Heavy Duty (K0824/K0825)
Weight Capacity	450 lbs.
Turning Radius³	22"
Overall Width³	25.4"
Overall Length³	40.5"
Maximum Speed²	Up to 4.88 mph
Ground Clearance³	3"
Front Casters	6"
Drive Wheels	13.75" perforated
Rear Casters	6"
Drivetrain	Two in-line HD motors
Braking System	Regenerative and electromechanical
Suspension Type	Full independent front and rear Active-Trac® ATX Suspension
High-back Seat Sizes	18" x 18-20" 20" x 18-20" 22" x 20-22" 24" x 22-24"
Seat-to-Floor Range	22-24"
Per Charge Range^{2,5}	Up to 13 miles at 450 lbs. Up to 21.5 miles at 200 lbs.
Standard Electronics	60A Dynamic Linx
Component Weights⁶	Base: 125.7 lbs. Seat: 57.8 lbs. Batteries: (2) 36 lbs. each;
Total Weight⁴	255.5 lbs.
Battery Requirements	(2) NF-22
Battery Charger	5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 13-month limited Electronics: 13-month limited Battery: 13-month limited

- ¹ Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain, battery type, battery condition and tire pressure.
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ Total weight includes base, seat and batteries.
- ⁵ Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.
- ⁶ Battery weight may vary based on manufacturer.



**PRIDE®
MOBILITY**

800-800-8586 • pridemobility.com



HCPSC codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

**Jazzy®
EVO 614HD**

Our most
advanced **heavy duty**
power chair yet!



800-800-8586 • pridemobility.com



The Jazzy EVO 614HD is an FDA Class II Medical Device.

THIS IS JAZZY EVOLVED

Height, width, depth
and angle adjustable
armrests

Removable
color-through
shrouds

Access to controller and
batteries without
the user needing to
exit the chair

Optional
swing-away
joystick

Depth-adjustable,
limited recline,
memory foam
high-back seat
with swivel

Heavy duty
seat post

13.75" perforated
drive wheels

FDA Class II
Medical Device

Front and rear independent
Active-Trac® ATX Suspension



High-back memory
foam seat



Large footplate with removable
rubber cover for easy cleaning



Front and rear independent
Active-Trac® ATX Suspension

FEATURES

- 450-lb. weight capacity
- 13.75" perforated mid-wheel drive tires
- Front and rear independent Active-Trac® ATX Suspension
- Removable matte-finish color-through shrouds
- Depth-adjustable, limited recline, memory foam heavy duty single seat post allows the user to swivel (max. seat size of 24" x 24")
- Access to controller and batteries without the user needing to exit the chair

AVAILABLE COLORS (MATTE)



Robin's
Egg Blue



Iceberg
Blue



Sugar
Plum



Black



White



Red

ACCESSORIES

- Attendant control bracket
- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag kit
- Lap belt
- Military patches
- Oxygen tank holder
- Rear basket
- Saddlebag
- Walker holder
- Weather cover
- 2.1A XLR USB charger
- Stealth® cushions: Zen/Essence/Glacial

OPTIONS

- Elevated leg rests
- Solid seat pan with limited recline and swivel
- Heel loops
- Swing-away joystick
- Swing-away leg rests
- Unoccupied transit
- Residual limb support