

Jazzy® EVO 614¹	Group 2: Standard (K0822/K0823)
Weight Capacity	300 lbs.
Turning Radius³	22"
Overall Width³	25.4"
Overall Length³	40.5"
Maximum Speed²	Up to 4.4 mph
Ground Clearance³	3"
Front Casters	6"
Drive Wheels	13.75" perforated
Rear Casters	6"
Drivetrain	Two in-line HD motors
Braking System	Regenerative and electromechanical
Suspension Type	Full independent front and rear Active-Trac® ATX Suspension
High-back Seat Sizes	16"x16-18" 18"x18-20" 20"x18-20"
Seat-to-Floor Range	22-24"
Per Charge Range^{2,5}	Up to 19 miles at 300 lbs. Up to 21.5 miles at 200 lbs.
Electronics	60A Dynamic Linx
Component Weights⁶	Base: 125.7 lbs. Seat: 57.8 lbs. Batteries: (2) 36 lbs. each;
Total Weight⁴	255.5 lbs.
Battery Requirements	(2) NF-22
Battery Charger	5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 13-month limited Electronics: 13-month limited Battery: 13-month limited

- ¹ Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain, battery type, battery condition and tire pressure. This specification can be subject to a variance of +10%, -5%.
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ Total weight includes base, seat and batteries.
- ⁵ Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities.
- ⁶ Battery weight may vary based on manufacturer.

**PRIDE®
MOBILITY**

800-800-8586 • pridemobility.com



HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

060926

**Jazzy®
EVO 614**

Always insist on **lithium!**



800-800-8586 • pridemobility.com



THIS IS JAZZY EVOLVED



FEATURES

- 3" ground clearance
- 13.75" perforated mid-wheel drive tires
- Front and rear independent Active-Trac® ATX Suspension
- Removable matte-finish color-through shrouds
- Depth-adjustable, limited recline, memory foam heavy duty single seat post allows the user to swivel
- Access to controller and batteries without the user needing to exit the chair

AVAILABLE COLORS (MATTE)



ACCESSORIES

- Attendant control bracket
- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag kit
- Lap belt
- Military patches
- Oxygen tank holder
- Rear basket
- Saddlebag
- Walker holder
- Weather cover
- 2.1A XLR USB charger
- Stealth® cushions: Zen/Essence/Glacial

OPTIONS

- Depth-adjustable solid seat pan, limited recline, high-back with swivel
- Elevated leg rests
- Swing-away leg rests
- Unoccupied transit
- Heel loops
- Residual limb support



High-back memory foam seat



Large footplate with removable rubber cover for easy cleaning



Front and rear independent Active-Trac® ATX Suspension