

Jazzy Select 2S-SS	Group 2: Standard (K0822)
Jazzy Select 2S-C	Group 2: Standard (K0823)
Weight Capacity	300 lbs.
Turning Radius	22.5"
Overall Width³	23.75"
Overall Length³	32.99"/44" (with footplate)
Maximum Speed²	Up to 4.3 mph
Ground Clearance³	.75" at motor bracket
Front Wheels	5" solid
Drive Wheels	10" solid
Rear Wheels	6" solid
Drivetrain	Two motor, mid-wheel drive®
Braking System	Regenerative and electromechanical
Suspension Type	Active-Trac® (ATS)
High-back Seat Sizes	16"x16-18", 18"x18-20", 20"x18-20"
Seat-to-Floor Range	21"-22"
Per Charge Range^{3,5}	Up to 19.4 miles
Standard Electronics	40A, PG GC3
Component Weights	Base 94 lbs. (heaviest piece) Seat 50.2 lbs. Batteries 25.2 lbs. each
Batteries^{4,6}	(2) 12 volt, U-1 batteries
Battery Charger	3.5 amp, off-board
Battery HCPCS	E2365
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

- 1 Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 AGM or gel cell type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

061226



PRIDE® MOBILITY

800-800-8586 • pridemobility.com



HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

Jazzy® SELECT

ACCESS YOUR WORLD WITH EASE



800-800-8586 • pridemobility.com



The Jazzy® Select is an FDA Class II Medical Device.

Jazzy® SELECT

Depth-adjustable,
limited recline,
fold-flat high
back seat

Intuitive controller is
easy to operate

Height, width, and angle
adjustable armrests

Jazzy® Red
color-through
shroud

Height and depth
adjustable foot
platform

Easy-access
free-wheel lever

FDA Class II Medical Device

FEATURES

- Tight turning radius of only 22.5"
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Articulating rear caster beam works in conjunction with the ATX Suspension for a smooth ride over uneven surfaces
- Easy-access free-wheel lever
- Color-through shroud keeps the shroud looking new by reducing the visual scrapes and chips

INCLUDED:

Deluxe contour, limited recline, high back with solid seat pan, available sizes:

- **16-20" (W) x 16-20" (D)**

Limited recline high back, available sizes:

- **16" (W) x 16" (D)**

OPTIONS

- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away footrests
- Swing-away joystick

ACCESSORIES

- Attendant control bracket
- Cane/crutch holder
- Cup holder
- Flag
- Lap belt (50", 70", 80", 90")
- Oxygen holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover (medium)
- XLR USB
- Stealth Products® cushions: Zen, Essence and Glacial



Easy battery access
from the side

Active-Trac® Suspension

PG GC3 Controller